




LAPAROSCOPY

INSUFFLATION LINE




ADVANTAGES

- Insufflation lines are designed to **withstand high gas flow rates** to allow **rapid, controlled** and **safe** insufflation
- Tubing is **flexible, non-obstructive, lightweight** and **highly resistant to plication**
- Device packaged in **double peelable sterilisation bag** and **traceability label**

REY MEDICAL SA: YOUR PARTNER IN SURGERY

• **QUALITY** • **SECURITY** • **ECONOMY**

Visit our website: www.reymedical.ch or call us: 022 743 24 10




Reliability

Extra-flexible tubing for anti-plication.


The polycarbonate filter and membrane with high hydrophobic, antibacterial and antiviral characteristics provides an excellent barrier to patient/insufflator cross contamination.

Filtration efficiency:
Bacteria: 99.999%
Particles and viruses: 99.99%



Compatibility

The various fittings ensure compatibility with all the insufflators. The wide range of fittings makes it possible to homogenise the equipment. The secure luer-lock fitting fits on all available trocars.




Ease of use

The ultra-light tubing avoids placing too much weight on the trocar or pneumo-peritoneum needle fitting.

Length of the tubing is 3m

Codes:



Connector 15 and 22 mm

LIF901522



Connector Push

LIF906220



Connector Toroidal

LIF906210



Connector Luer Lock

LIF906250