




LAGIS[®]
MINIMAL INVASIVE INNOVATED

LapShield[™]
Wound Retractor/Protector



LapShield™

Wound Retractor

Elastic Abdominal Retractor for Minilaparotomy

- Maximize exposure, minimize incision size
- Reduce the risk of wound infections
- Offer ultimate versatility



360° Atraumatic Retraction

- Distributes force evenly for optimal field of view
- Allows for maximum exposure with a minimum incision size
- Soft materials and evenly dispersed tension promote patient post-operative comfort

360° Protection


- Significantly decreases wound infection
- Allows visualization of wound margins
- Maintains moisture at the incision site
- Helps protect wound from contamination while helping to keep wound edges moist

Ultimate Versatility

- Unique design facilitates rapid and effortless set-up
- Multiple sizes accommodate a wide range of surgeries

Easily Placed and Removed

- Inserted in seconds
- Flexible ring slips gently into incision
- Removed with one finger



www.lagis.com.tw

©2013 LAGIS-L09-C-2018.01

Specifications subject to change without notice.

Product Code	Incision Length (mm)	Quantity
WR-30US	10~20	6 /box
WR-60ES	20~40	6 /box
WR-80S	30~60	6 /box
WR-110M	50~90	6 /box
WR-160L	90~140	6 /box
WR-190EL	110~170	6 /box

Designed for soft tissue retraction, Lagis LapShield™ provides unparalleled exposure without the trauma associated with prolonged point retraction.

service@lagis.com.tw

Not made with natural rubber latex.

REY MEDICAL SA: YOUR PARTNER IN SURGERY

• **QUALITY** • **SECURITY** • **ECONOMY**

Visit our website: www.reymedical.ch or call us: 022 743 24 10



LapBase™

Advanced Single Access Platforms



**Reduced Port Endoscopic Surgery
with LapBase™ to create Better Performance**

LAGIS®

LapBase™ Cap

Manipulation enhanced by freely trocar positioning

- Offer an abdominal-like platform for single incision procedures
- Able to manipulate standard laparoscopic instrument
- Stable platform and maintain pneumoperitoneum with every maneuver
- Direct transparent visualization



LapShield™ Wound Protector

Reduced Port Access Platforms

- Provide 360 degrees protection
- Allows extracorporeal anastomosis and specimen retrieval
- Maximum exposure with minimum incision size
- Distribute force evenly for optimal field of view



LapShield™ Wound Protector

Product code	Incision Length (mm)
WR-60ES	20 - 40
WR-80S	30 - 60
WR-110M	50 - 90

Packaging: 6EA/ box

LapBase™ Cap

Product Code	Description	Incision Length (mm)
CP-60ES	Compatible with WR-60ES	20 - 40
CP-80S	Compatible with WR-80S	30 - 60
CP-110M-A	Compatible with WR-110M	50 - 90

Packaging: 6EA/box



LapBase™ Kit

Product Code	Description
LB-6001	WR-60ES + CP-60ES + TRO-12V1 (12mm ClearView™ Optical Trocar, 108mm length)
LB-6002	WR-60ES + CP-60ES + TRO-12C1 (12mm Dilating Tip Trocar, 108mm length)
LB-8001	WR-80S + CP-80S + TRO-12V1 (12mm ClearView™ Optical Trocar, 108mm length)
LB-8002	WR-80S + CP-80S + TRO-12C1 (12mm Dilating Tip Trocar, 108mm length)
LB-1100	WR-110M-A + CP-110M
LB-1101	WR-110M-A + CP-110M + TRO-12V1 (12mm ClearView™ Optical Trocar, 108mm length)
LB-1102	WR-110M-A + CP-110M + TRO-12C1 (12mm Dilating Tip Trocar, 108mm length)

Packaging: 1EA/box

Caution: Always use a needle like conical tip to create the first insertion within the boundary line, do not penetrate the LapBase Cap with trocar directly.

REY MEDICAL SA: YOUR PARTNER IN SURGERY

• QUALITY • SECURITY • ECONOMY

Visit our website: www.reymedical.ch or call us: 022 743 24 10



- Offer an abdominal-like platform for single incision procedures
- Able to manipulate standard laparoscopic instrument
- Stable platform and maintain pneumoperitoneum with every maneuver
- Direct transparent visualization LapShield



Product Code	Compatible Wound Retractor	Incision Length (mm)
CP-60ES	WR-60ES	20 – 40
CP-80S	WR-80S	30 – 60
CP-110M	WR-110M-A	50 – 90

Department	Single Port Procedure	Laparoscopic	Robotic
Bariatric	Gastric Bypass	CP-80S	-
	Sleeve Gastrectomy	CP-60ES	CP-110M
Colon & Rectal	Colectomy	CP-60ES	-
	Right Hemicolectomy	CP-80S	CP-80S
General Surgery	Appendectomy	CP-60ES	-
	Cholecystectomy	CP-60ES/CP-80S	CP-60ES/CP-80S
	Inguinal Hernia Repair	CP-60ES/CP-80S	-
	Whipple	-	CP-110M
	Splenectomy	CP-60ES	-
	Distal Pancreatectomy	CP-80S	CP-110M
	Spleen-Preserving Distal Pancreatectomy	CP-60ES/CP-80S	CP-80S
OB/GYN	Ovarian Cystectomy	CP-60ES/CP-80S	CP-110M
	Total Laparoscopic Hysterectomy(TLH)	CP-60ES	CP-80S
	Laparoscopy-Assisted Vaginal Hysterectomy	CP-60ES/ CP-80S	-
	Myomectomy	CP-60ES	CP-80S
Urology	Total Extraperitoneal (TEP) Herniorrhaphy	CP-80S	-
	Radical Nephrectomy	CP-80S	CP-80S/ CP-110M
	Prostatectomy	CP-80S	-

REY MEDICAL SA: YOUR PARTNER IN SURGERY

• **QUALITY** • **SECURITY** • **ECONOMY**

Visit our website: www.reymedical.ch or call us: 022 743 24 10

LapShieldTM Wound Retractor

Department	Procedure	WR-30US 1~2cm incisions	WR-60ES 2~4cm incisions	WR-80S 3~6cm incisions	WR-110M 5~8cm incisions	WR-160L 9~14cm incisions	WR-190EL 11~17cm incisions
BARIATRIC	Lap Gastric Bypass		v	v			
BARIATRIC	Lap Sleeve Gastrectomy	v					
BARIATRIC	Open Gastric Bypass					v	v
COLON & RECTAL	Open Colectomy					v	v
COLON & RECTAL	Lap Colectomy			v	v		
COLON & RECTAL	Stoma Formation	v					
General Surgery	Appendectomy			v	v		
General Surgery	Inguinal Hernia Repair		v	v			
General Surgery	Pancreatectomy					v	v
General Surgery	Splenectomy					v	v
General Surgery	Thyroidectomy		v	v			
OB/GYN	Bilateral Salpingo Oophorectomy		v	v			
OB/GYN	Cesarean Section					v	v
OB/GYN	Lap Hysterectomy			v	v		
OB/GYN	Mini-Laparotomy			v	v		
OB/GYN	Myomectomy			v	v		
OB/GYN	Postpartum Tubal Ligation	v	v				
OB/GYN	Total Abdominal Hysterectomy			v	v	v	
THORACIC	Mitral Valve Repair/Replacement			v	v		
THORACIC	Thoracotomy			v	v		
THORACIC	Video-Assisted Thoracoscopic Surgery (VATS)	v	v	v			